



Y-90 SPECT maximum likelihood image reconstruction with a new model for tissue-dependent bremsstrahlung production

- 2017 SNM Annual Meeting -


Hongki Lim, Yuni Dewaraja, Scott Wilderman, Jeffrey Fessler
University of Michigan, Ann Arbor

Disclosure

- Yuni Dewaraja and Jeffrey Fessler are consultants for MIM software Inc., Cleveland, Ohio

Objectives


- Material dependence of the bremsstrahlung generation
 - Bremsstrahlung production is **proportional to Z^2** of absorber



Bremsstrahlung generation illustration


Objectives

- Y-90 SPECT
 - Existing reconstruction methods do not account for tissue-dependency




True image with equal activity conc. in bone & marrow

Reconstruction
→




Reconstructed image



CT

Methods

- Monte Carlo simulation
 - Bremsstrahlung generation in bone is ~2 times higher than in other tissues



- Simulation Software: EGS5 [1]
- Simulated for Bone, Tissue, Marrow, Lung in an infinite media

➡ Same amount of activity in different tissue will produce different numbers of bremsstrahlung photons

Methods

- New system matrix
 - Incorporated the tissue-dependent bremsstrahlung generation probability (based on CT) into the reconstruction system matrix :

$$y \approx \text{Poisson}(Ax + s)$$

Where, $A = A_{\text{attn+psf}} B_{\text{prob}}$

$$B_{\text{prob}} = \text{diag}(b)$$

$$b_j = q_{\text{tissue}} * (1 - \text{BVF}_j) + q_{\text{bone}} * (\text{BVF}_j)$$

$$q_{\text{bone}} \simeq 2 * q_{\text{tissue}}$$

BVF_j : Bone-Volume Fraction of j th voxel in SPECT

- Modification to standard system matrix:

$$A = \begin{bmatrix} a_{11} & a_{12} & \dots & a_{1n} \\ a_{21} & a_{22} & \dots & a_{2n} \\ \dots & \dots & \dots & \dots \\ a_{n1} & a_{n2} & \dots & a_{nn} \end{bmatrix} \begin{bmatrix} b_1 & 0 & \dots & 0 \\ 0 & b_2 & \dots & 0 \\ \dots & \dots & \dots & 0 \\ 0 & 0 & \dots & b_n \end{bmatrix}$$

Standard system matrix

Our modification

Methods


- BVF example:

B	B	B	B	B
B	B	B	B	B
T	B	B	B	T
T	B	B	B	T
B	B	B	T	T

25 voxels in CT image

* B: Bone, T: Tissue

→ Determined by thresholding




1 voxel in SPECT image

Results

- Measurement setup
 - Same Y-90 concentration in syringes with water and bone equivalent solution


	Water Syringe	Bone Syringe
Injected Y-90 Amount	6.94mCi	6.81mCi
Solution	Water	Bone equivalent solution [2] (K ₂ HPO ₄ , Dipotassium Hydrogen Phosphate)
Solution Density	~ 1.0 g/cm ³	~ 1.68 g/cm ³
Solution Amount	40 cc	40 cc

- CT Image:




Results

- Bremsstrahlung yields from our EGS5 simulations
- Attenuation coefficients from [2]



Photon generation vs Photon Energy




Attenuation coefficient vs Photon Energy

Results


- Commercial PET and SPECT OSEM

PET
(consistent w/ TRUE)



- Our in-house SPECT OSEM reconstruction

Standard method




New method
(consistent w/ TRUE)



Results


- Profile comparison:




- Slice comparison



True Image




Standard Reconstruction




New Reconstruction

Results


- Simulation setup
 - XCAT phantom slices from spine to upper femur
 - Attenuation map at three slices:



Spine




Pelvis



Upper femur


- Set amount of activity in bone and marrow region as 1:1, 3:1, 1:3




True Activity Map (1:1 case) at Pelvis

Results


- Simulation result (ratio = 1:1)
 - Qualitative result



True Image



New Reconstruction




Standard Reconstruction


Quantitative Results

- Simulation result (ratio = 1:1)


Count Recovery @ Bone




Count Recovery @ Marrow




— New model — Standard model — TRUE — New model — Standard model — TRUE



[Lumbar]



[Pelvis]



[Upper femur]

Quantitative Results

- Simulation result
 - When activity ratio between bone:marrow is 3:1 and 1:3

Uptake ratio (Bone:marrow)	Model	ROI	Count Recovery @ Marrow	Count Recovery @ Bone
3:1	New model	Average	169.4	69.4
	Standard model		217.9	110.6
1:3	New model		81.6	109.6
	Standard model		92.6	160.3

Summary / Conclusion

- Demonstrated the tissue-effect in simulation and measurement
- Demonstrated proof of concept of reconstruction incorporating tissue-dependent bremsstrahlung generation
- Potentially applicable to Y-90 therapies such as radio immunotherapy and synovectomy

Acknowledgements

- NIH R01EB022075
- MIM for providing research license
- UM Radiology:
Allen Brooks, David Hubers, Jeremy Niedbala

Thank you.

Contact: hongki@umich.edu

