

Joint Dual Photopeak Image Reconstruction in Lu-177 SPECT

2018 EANM

Hongki Lim, Jeffrey Fessler, Yuni Dewaraja
University of Michigan, Ann Arbor



EANM Disclosure Statement

1) I or one of my co-authors hold a position as an employee, consultant, assessor or advisor for a pharmaceutical, device or biotechnology company. **If yes, please specify name/position/company:**

Nothing to disclose

2) I or one of my co-authors receive support from a pharmaceutical, device or biotechnology company. **If yes, please specify name/position/company/which project and whether support is in kind or monetary:**

Yuni Dewaraja is paid consultants for MIM Software Inc., Cleveland, USA

3) I or one of my co-authors hold property rights/patents for (radio)pharmaceuticals, medical devices or medical consulting firms. **If yes, please specify name/position/company:**

Nothing to disclose

4) I or one of my co-authors have written articles for (radio)pharmaceutical, medical device, biotechnology or consulting companies during the last 5 years. **If yes, please specify name/position/company/article/ journal and co-authors:**


Nothing to disclose



Energy Spectrum of ^{177}Lu

Simulation with point source in water phantom & ME collimator

- General recommendation¹
 - ME collimator + 208 keV
- Another option
 - ME collimator + 113 & 208 keV



1. Ljungberg et al. "MIRD pamphlet no. 26" J Nucl Med 57.1 (2016): 151-62.



Energy Spectrum of Lu-177


General recommendation

- ME collimator + 208 keV
→ Clinically adopted option

Another option

- ME collimator + 113 & 208 keV

Simulation with point source in water phantom & ME collimator



Energy Spectrum of Lu-177

General recommendation


- ME collimator + 208 keV

Another option

- ME collimator + 113 & 208 keV

→ Good to use when # of counts is insufficient (imaging at later time point or pretherapy diagnostic imaging)


Simulation with point source in water phantom & ME collimator



Energy Spectrum of Lu-177

Simulation with point source in water phantom & ME collimator

- Compared to 208-keV, 113-keV has
 - ~60% of Yield
 - ~70% Primary Counts (cps/MBq) for ME collimator



Options for using both photopeaks - 113 & 208 keV

Two options


1) Reconstruct two images **separately** and make sum

$$\begin{aligned}y_1 &\sim \text{Poisson}(A_1 x_1 + s_1), \hat{x}_1 = \operatorname{argmin}_{x_1} f(x_1; y_1, s_1) \\y_2 &\sim \text{Poisson}(A_2 x_2 + s_2), \hat{x}_2 = \operatorname{argmin}_{x_2} f(x_2; y_2, s_2) \\&\hat{x} = \hat{x}_1 + \hat{x}_2\end{aligned}$$

2) Reconstruct one image **jointly**

$$\begin{aligned}y &= \begin{bmatrix} y_1 \\ y_2 \end{bmatrix} \sim \text{Poisson}\left(\begin{bmatrix} A_1 \\ A_2 \end{bmatrix} x + \begin{bmatrix} s_1 \\ s_2 \end{bmatrix}\right) \\&\hat{x} = \operatorname{argmin}_x f(x; y, s)\end{aligned}$$

- $y_1(y_2)$ is main energy window and $s_1(s_2)$ is scatter energy window acquisition around 208(113)-keV
- $f(x)$ denotes negative Poisson likelihood function
- For scatter estimation, we used triple-energy-window (TEW) technique



Options for using both photopeaks - 113 & 208 keV

Two options

1) Reconstruct two images **separately** and make sum


$$\begin{aligned}y_1 &\sim \text{Poisson}(A_1 x_1 + s_1), \hat{x}_1 = \operatorname{argmin}_{x_1} f(x_1; y_1, s_1) \\y_2 &\sim \text{Poisson}(A_2 x_2 + s_2), \hat{x}_2 = \operatorname{argmin}_{x_2} f(x_2; y_2, s_2) \\&\hat{x} = \hat{x}_1 + \hat{x}_2\end{aligned}$$

2) Reconstruct one image **jointly**

$$\begin{aligned}y &= \begin{bmatrix} y_1 \\ y_2 \end{bmatrix} \sim \text{Poisson} \left(\begin{bmatrix} A_1 \\ A_2 \end{bmatrix} x + \begin{bmatrix} s_1 \\ s_2 \end{bmatrix} \right) \\&\hat{x} = \operatorname{argmin}_x f(x; y, s)\end{aligned}$$

→ More aligned with actual Poisson nature of multi-window acquisition

- $y_1(y_2)$ is main energy window and $s_1(s_2)$ is scatter energy window acquisition around 208(113)-keV
- $f(x)$ denotes negative Poisson likelihood function
- For scatter estimation, we used triple-energy-window (TEW) technique



What factors need to be considered for joint reconstruction?

□ System Model

$$\begin{bmatrix} y_1 \\ y_2 \end{bmatrix} \sim \text{Poisson} \left(\begin{bmatrix} A_1 \\ A_2 \end{bmatrix} x + \begin{bmatrix} S_1 \\ S_2 \end{bmatrix} \right),$$

where

$$a_{ijk} = \alpha_k p_{ijk} g_{ijk}$$

a_{ijk} : element of system model (A_k) at i -th row and j -th column in k -th energy window

α_k : Relative camera **detection efficiency** \times branching ratio

p_{ijk} : depth-dependent and energy-dependent effect of **attenuation**

g_{ijk} : depth-dependent and energy-dependent **3D-collimator detector response**



Experimental setting


□ 2 measurements + 1 simulation + 1 patient scan

- I. Sphere phantom measurement (4 scans)
- II. Liver-torso phantom measurement (4 scans with high/low-counts)
- III. XCAT phantom simulation (10 realizations with high/low-counts)
- IV. Patient scan




Physical phantom measurement 1


- ☐ Sphere phantom measurement (4 scans to assess noise)
 - Sphere-to-background ratio was **6:1** , 30 minutes acquisition
 - Main window width: 20%, scatter window width: 10 %




208-keV main window




208-keV scatter window



113-keV main window



113-keV scatter window



Physical phantom measurement 2


□ Liver-torso phantom measurement

- 4 scans to assess noise:

- 4 scan of 35 min (mimicking Day3 post ^{177}Lu -DOTATATE)

- 4 scans of 3 min (mimicking Day3 pre-therapy diagnostic imaging)


- Lesion-to-liver ratio was 4:1




XCAT simulation

- XCAT coupled with SIMIND² simulation (10 realizations to assess noise)
 - Realistic activity distribution for ¹⁷⁷Lu-DOTATATE


Organ	Relative activity concentration
Lesion in liver	1.00
Spleen	0.59
Kidney	0.48
Liver	0.16
Gall bladder	0.02
Renal pelvis	0.02




208-keV main




208-keV scatter




208-keV main




208-keV scatter




113-keV main



113-keV scatter



113-keV main



113-keV scatter

Post-therapy Day1 count-level

Pre-therapy Day1 count-level




2. Ljungberg M. Monte Carlo Calculation in Nuclear Medicine: Application in Diagnostic Imaging, 2nd edn, 2012.

^{177}Lu -DOTATATE Patient scan


Day5 acquisition

- 1 clinical scan with 208-keV (Scanned for 25 minutes)
- 3 research scans with both photopeaks (Scanned for 2 minutes to mimic pre-therapy)

208-keV main window




208-keV scatter window




25 minutes scan at Day 5


208-keV main window




208-keV scatter window



113-keV main window



113-keV scatter window



2 minutes scan at Day 5



Reconstruction and Evaluation methods

- Reconstruction details (for both single and joint)
 - Applied Gaussian smoothing to TEW scatter estimate
 - EM algorithm + 50 iteration (1 subset)

□ Evaluation metrics

I. Activity-recovery


$$\text{Activity recovery (\%)} = \frac{\text{Estimated mean activity in VOI}}{\text{True mean activity in VOI}} \times 100 \%$$

II. Contrast-recovery




$$\text{Contrast recovery (\%)} = \frac{\frac{\text{Estimated mean counts in VOI}}{\text{Estimated mean counts in BKG}} - 1}{\text{True activity concentration ratio between VOI and BKG} - 1} \times 100 \%$$

III. Image-ensemble-noise (variability across multiple acquisitions)

$$\text{Image ensemble noise (\%)} = \frac{\sqrt{\frac{1}{J_{\text{BKG}}} \sum_{j \in \text{BKG}} \left(\frac{1}{M-1} \sum_{m=1}^M (\hat{x}_m[j] - \frac{1}{M} \sum_{m'=1}^M \hat{x}_{m'}[j])^2 \right)}}{\frac{1}{J_{\text{BKG}}} \sum_{j \in \text{BKG}} \frac{1}{M} \sum_{m=1}^M \hat{x}_m[j]} \times 100 \%$$



Calibration factor for absolute quantification





1. NIST recommendation³ (0.9 % uncertainty)
 - 3 mL ¹⁷⁷Lu in a 10 mL Schott vial: CRC-15R dial setting 449x10
2. Transferring Calibration to a new geometry (10 mL syringe)
 - Activity in syringe = activity in vial - residual
3. With the syringe in the dose calibrator adjust setting
 - Setting to get correct reading for 3 mL in syringe: 480 x 10

- Imaging for calibration factor
 - planar with 20% main window & 10% TEW SC
- 208 keV (head1): 12.9 cps/MBq
(head2): 13.4 cps/MBq
→ (within 1 % of manuf. specified value for 5/8" xtal)
- 113 keV (head1): 9.2 cps/MBq (70% of above)
(head2): 9.3 cps/MBq

3. Bergeron DE, Cessna JT. Nucl Med Commun. 2018 Jun;39(6):500-504






Result: Sphere phantom measurement



Profile across cold spot




Result: Sphere phantom measurement




Result: Liver-torso phantom measurement (axial view)

High Counts →


Single 208-keV



Single 113-keV




Joint




Low Counts →



Single 208-keV



Single 113-keV



Joint



Result: Liver-torso phantom measurement

35 minutes scan


	Noise	Contrast Recovery		CRNR	
		29ml ovoid	16ml sphere	29ml ovoid	16ml sphere
208-keV	11.4	75.3	70.4	6.6	6.2
113-keV	16.4	72.4	66.6	4.4	4.1
Joint	9.5	74.4	69.1	7.8	7.2

3 minutes scan

	Noise	Contrast Recovery		CRNR	
		29ml ovoid	16ml sphere	29ml ovoid	16ml sphere
208-keV	36.4	65.2	63.7	1.8	1.8
113-keV	48.5	58.7	51.8	1.2	1.1
Joint	29.3	62.1	58.1	2.1	2.0




Result: XCAT phantom simulation



Post-therapy Day1 count-level

	Noise	Contrast Recovery		CRNR	
		29ml ovoid	16ml sphere	29ml ovoid	16ml sphere
208 keV	10.1	76.3	92.9	7.5	9.2
113 keV	12.9	73.3	90.4	5.7	7.0
Joint	8.0	75.0	91.8	9.4	11.5




Pre-therapy count-level

	Noise	Contrast Recovery		CRNR	
		29ml ovoid	16ml sphere	29ml ovoid	16ml sphere
208 keV	65.2	75.8	91.7	1.2	1.4
113 keV	85.5	68.7	91.5	0.8	1.1
Joint	52.2	72.8	91.6	1.4	1.8



Result: Patient scan



Recon (unit: Bq/ml)	RMSD*	Noise*
208-keV	1.47e4	332.9
113-keV	1.90e4	278.6
Joint	1.45e4	260.4

* RMSD relative to long scan

* Image ensemble noise between 3 repeat short scans



Summary

- Joint reconstruction w/ E-dependent system model & TEW SC is easy to implement
- Compared to single photopeak (208-keV) reconstruction, joint reconstruction
 - Gave comparable activity/contrast recovery
 - Improved noise significantly even for high-counts data
 - Important for voxel-level dosimetry
 - Improved contrast recovery to noise ratio: 16.7-27.3%
- Joint reconstruction can be beneficial for reducing scan time

Improvement (%) of Contrast Recovery to Noise Ratio over single 208-keV

		Experiments				
		Sphere phantom	Liver phantom	XCAT		
Count-level		High	Low	High	Low	
Improvement(%)		20.9-23.2	16.7-17.5	17.3-21.5	25.4-25.8	24.7-27.3



We acknowledge

- U of Michigan Radiology Seed Grant
- Jeremy Niedbala, Gerrid Rosebush
for all help in measurements and calibration

Thank You

Contact: hongki@umich.edu

